

Action



Pressure relieving cushion and cover compliance: Monthly Audit Tool for Care Homes

This Audit Tool aims to identify pressure relieving cushions and covers in poor condition. Cushions and covers in a poor condition are known to be a source of healthcare associated infection.

- Each pressure relieving cushion and cover should be reviewed on a regular basis, which should be based upon a local risk assessment, e.g. monthly. Complete sections 1 and 2 below for each pressure relieving cushion cover and cushion.
- Cushions should be enclosed in a waterproof cover, preferably with an integral zip fastener, to facilitate inspections of the surfaces.
- For audit purposes, all cushion covers should be numbered for identification, rather than identifying the cushion by the location, as cushions are often moved between locations.
- In the event of cushion or cover failure, action plans must be drawn up locally, reviewed and monitored.
- Completed audit tools should be kept locally for good practice assurance and as evidence for CQC inspections.

Name of establishment	Audit completed by	
	Job title	
Unit name (if applicable)	Date	

	1. Criteria for pressure relieving cushion covers								
Cushion number and location	Is there a breach in the integrity of the cushion cover, e.g. torn or damaged?		Is the cushion cover fastening compromised, e.g. is the zip or any other cover fastening device broken?		Does the cushion cover have any staining outside or inside that cleaning cannot remove?		If the answers to any of the questions are yes, the pressure relieving cushion cover has failed and must be replaced		
	Yes	No	Yes	No	N/A	Yes	No	Fail	Pass
1.									
2.									
3.									
4.									
5.									
6.									
7.									
8.									
9.									
10.									

Name of establishment		Audit completed by	
Name of establishment		Job title	
Unit name (if applicable)		Date	

	2. Criteria for pressure relieving cushions							
Room number/cushion location	Undo the removable co on either side? (Not applicable to seale			Does the cushion hav odour?	e an offensive	If the answers to any of the questions are yes, the pressure relieving cushion has failed and must be replaced		
	Yes	No	N/A	Yes	No	Fail	Pass	
1.								
2.								
3.								
4.								
5.								
6.								
7.								
8.								
9.								
10.								

© Harrogate and District NHS Foundation Trust, Community Infection Prevention and Control <u>www.infectionpreventioncontrol.co.uk</u>