



Mattress and mattress cover compliance: Audit Tool for Care Homes

This Audit Tool aims to identify mattresses and mattress covers in poor condition. Mattresses and covers in a poor condition are known to be a source of healthcare associated infection. This 'Mattress and mattress cover Audit Tool' has been produced to identify damaged covers and associated interior mattress contamination. A frequent inspection regime should be established for all mattresses. Refer to 'Safe management of care equipment Policy for Care Home settings'.

Action

- Each mattress and mattress cover should be reviewed on a regular basis, which should be based upon a local risk assessment, e.g. weekly, monthly, and whenever a room is vacated. Following the instructions below, complete sections 1 and 2 for each mattress and mattress cover.
- Mattresses should be enclosed in a waterproof cover, preferably with an integral zip fastener, to facilitate inspections of the surfaces.
- All mattress covers should be numbered for identification, rather than identifying the mattress by the room number for audit purposes, as mattresses are often moved between rooms. It may not be possible to label some mattress covers, i.e. special mattresses. An appropriate system must be in place to identify these mattresses.
- In the event of mattress or mattress cover failure, an 'IPC Action Plan' should be produced and reviewed regularly.
- Completed Audit Tools and Action Plans should be kept locally for good practice assurance and as evidence for CQC inspections.

Instructions for checking a mattress and cover

Hand compression assessment

- 1. Where appropriate, adjust the height of the bed so that it is at the same level as the tester's hip.
- 2. Link hands to form a fist and place them on the mattress.
- 3. Keep elbows straight and lean forward, applying the full body weight to the mattress.
- 4. Repeat the hand compression at intervals along the entire length of the mattress.
- 5. The mattress should be condemned if the base of the bed can be felt at any point.

Water penetration test (only for sealed mattress covers without a zip)

Water penetration testing should only be carried out on mattresses that have a sealed cover, i.e. no zip, which make it impossible to inspect the core mattress. This is intended as a quick positive check to determine whether or not the mattress cover is permeable to fluids.

- 1. Using the fist, press the mattress cover to form a shallow well in the centre of the mattress where it would have high usage or on any areas of concern.
- 2. Pour a small volume of tap water (approximately 75 ml) in the well and agitate the area for about one minute.
- 3. Visually inspect for loss of water/water penetration which indicates failure of the test.
- 4. Mop up the water.

Turning mattresses

Many new mattresses do not require turning. However, if the mattress does require turning, this should be carried out in accordance with the manufacturer's instructions.

Name of establishment Unit name (if applicable)			Audit completed by											
			Job title											
			Date											
		1. Criteria for mattress covers												
	Mattress number/ identifier	Is there a breach in the integrity of		Removable mattress covers: is the mattress cover fastening			Non-removable mattress covers, e.g. no zip: did the			Does the mattress cover have any		If the answers to any of the		
Room number		the mattre	ess cover,	compromised, e.g. is the zip or			cover fail the 'Water			staining outside or		mattress cover questions are yes, the cover has		
		e.g. torn or damaged?		any other cover fastening device broken?			penetration test'? (see page 1)			inside that cleaning cannot remove?		failed and must be replaced		
		Yes	No	Yes	No	N/A	Yes	No	N/A	Yes	No	Fail	Pass	
1.														
2.														
3.														
4.														
5.														
6.														
7.														
8.														
9.														
10.														
11.														
12.														
13.														
14.														
15.														
16.														
17.														
18.														
19.														
20.														

Name of establis					Audit co	ompleted by	,					
Name of establis							Job title					
Unit name (if applicable)			Date									
	_											
Room number	Mattress number/ identifier	Undo the remattress we either side?	a for mattre emovable cover, soiled or soiled	er. Is the tained on able to	Does the m have an offe odour?		compressio	etress fail the ' on assessmen ot applicable to)	ť? (see	If the answers to any of the questions are yes, the mattress has failed and must be replaced		
		Yes	No	N/A	Yes	No	Yes	No	N/A	Fail	Pass	
1.												
2.												
3.												
4.												
5.												
6.												
7.												
8.												
9.												
10.												
11.												
12.												
13.												
14.												
15.												
16.												
17.												
18.												
19.												
20.												

Available free to download, go to www.infectionpreventioncontrol.co.uk (© Harrogate and District NHS Foundation Trust, Community Infection Prevention and Control www.infectionpreventioncontrol.co.uk